

CARECLEAN CITRIC ACID

Concentrated Citric Acid powder for tank cleaning.

DESCRIPTION

Careclean Citric Acid is concentrated citric acid powder. Careclean Citric Acid can be used to remove hardness stains, passivate stainless steel surfaces and is an excellent alternative to strong acid cleaners.

ADVANTAGE

- Mild Acid
- Safe removal of Hardness stains
- ASTM method for passivation of stainless steel surfaces

APPLICATIONS

Tank Cleaning—Removal Hardness stains or tank brightener

Prepare a solution of 1% Careclean Citric Acid and recirculate for 3-4 hours at a temperature of 60 °C.

Passivation of stainless steel

Prepare a solution of 4-20% Careclean Citric Acid and recirculate at a temperature of 60 °C for at least 1 hour. Rinse with DI or low chloride water.

DOSAGE

Recirculation	1,0 %
Passivation	4,0 –20,0 %

APPROVALS

- IMO approved in accordance with MEPC 590 Guidelines
- For use on MarineLINE tank coatings, following restriction applies: dosing 0.1-4% at max. 35 °C

For detailed information on safety and health, please refer to the Material Safety Data Sheet MSDS and/or product label.

- Concentrated Citric Acid powder
- Mild pH tank brightener
- Removal of Hardness stains
- Passivation of Stainless steel ASTM A967-05
- Approvals:
 - > IMO



Revision	2.0	Document number	PDS-08-121
Date of issue	28-07-2023	Page	1 of 1

MARINE CARE BV
Oude Maasweg 35
Port number 4005
3197 KJ Rotterdam
The Netherlands
T. +31 (0)10 2950342
F. +31 (0)10 2950345
E. operations@marinecare.nl
W. www.marinecare.nl